



FITTING INSTRUCTIONS

EZ LACER

*Fits left or
right ankle*

1. Open both the front entry flap and the elastic foot control strap. Loosen laces on the front entry strap.
2. Insert foot into brace so that the heel rests in rear opening. The front end of the brace rests approximately midfoot.
3. Close both the front entry flap and the elastic foot control strap securely. Closures will provide some pressure but should not be uncomfortable or too snug.
4. Tighten and tie laces to provide a secure fit that envelops the lower leg.

*Soft tissue grip is key. For maximum control you may need to choose one size smaller brace.

For more information call 1-800-800-2896