



FITTING INSTRUCTIONS



Fits left or right ankle

- Volt™ will require a brief break-in period.
- For maximum support, wear a sneaker or a lace-up shoe.
- Wear a comfortable, absorbent sock at least calf-height.

- 1.** Place Volt™ inside your shoe with the logo label positioned so that it will be on the outside of your ankle.



- 2.** If your shoe has a removable insole or orthotic, place Volt™ under it for additional comfort. If doing so, we highly recommend removing the white heel pad.

- 3.** Place your foot inside the shoe. Make sure to adjust the pivot points to be in line with your ankle bones.



- 4.** Place the small posterior strap comfortably and without tension around the back of your ankle and secure it as low as possible to the Velcro® on the inside.

- 5.** Secure the large fastening straps firmly around the ankle and tie shoelaces. Now you are ready and all charged up!



For more information call 1-800-800-2896